

Part No.	8050H1500S
Description	Disposable LiquiThief – 1500mm
Material of Construction:	
Tip	HDPE (virgin, FDA approved – CFR 177.1520)
Outer Tube	HDPE (virgin, FDA approved – CFR 177.1520)
Inner Rod	HDPE (virgin, FDA approved – CFR 177.1520)
Clutch	HDPE (virgin, FDA approved – CFR 177.1520)
Piston Housing	HDPE (virgin, FDA approved – CFR 177.1520)
Piston Seal	Silicone Rubber (virgin, FDA approved – CFR 177.2600)
Handle	HDPE (virgin, FDA approved – CFR 177.1520)



Method of Construction:	Ultrasonic welding
Nominal Length	
Fully closed:	1495mm
Fully extended:	2830mm
Diameter of Outer Tube:	21mm
Diameter of hole in tip:	1.5mm
Nominal Weight of Sampler:	240g
Assembly & Packing Environment:	Class 100,000 Medical Cleanroom
Individually Bagged?	Yes (heat sealed PE bag)
Method of Sterilisation:	Electron Beam
Number of Samplers per Box:	20
Recommended Storage Conditions:	Dry and ambient temperature