Disposable Pod Powder Sampler Minimize Sampling Error



The Disposable Pod Sampler is the next generation in sampling technology. It combines a single use disposable Pod with a reusable stainless steel sample thief.

The sampler is designed to minimize sampling error, which can occur when a sample is transferred from one container to another. With the Accurate Pod the sample remains in the pod until it is needed in the laboratory.

Ideal for sampling free flowing powders and granules.




- Insert the Pods into the sample thief
- Take sample
- After sampling click the lids onto the Pods
- Take the Pods to the laboratory
- Slide off the lid to access the sample





It's as easy as 1, 2, 3



Pods

- Made from FDA, EU and EC compliant HDPE
- Initially available in 1ml and 5ml sizes. Other sizes to follow please contact us for details
- Supplied with lid (lid supplied loose in the bag)
- Individually bagged
- Made in a Class 8 (class 100,000) cleanroom
- Available pre sterilised (gamma irradiation)

Disposable Pod Sample Thief

Part No.	Nominal Volume	Sterile	Material	Pack Size
675 8360T10	1ml	No	HDPE	50
675 8360T10S	1ml	Yes	HDPE	50
675 8360T50	5ml	No	HDPE	50
675 8360T50S	5ml	Yes	HDPE	50

QUICK QUOTE

- Made from 316L stainless steel
- Available in a range of different lengths
- Available from stock
- 25mm Diameter

Single Port

Part No.	Sampler Length	No. of Ports
675 1360A601	600mm	1
675 1360A1001	1000mm	1
675 1360A1501	1500mm	1
675 1360A2001	2000mm	1

ONLINE PRICING

• Samplers with a bespoke number of ports available on request

Multi Port

316L		Stock	Download
	Video	Items	certificates

ONLINE PRICING

Part No.	Sampler Length	No. of Ports
675 1360A603	600mm	3
675 1360A1003	1000mm	3
675 1360A1503	1500mm	3
675 1360A2003	2000mm	3

QUICK QUOTE

QAQC LAB 589 Rappahannock Drive White Stone Va 22578 TEL (866) 244-1578 emails: sales@qclabequipment.com