

Multi Thief

Multilevel Powder Sampler with variable Sample Volumes

The Multi Thief is ideal for taking multilevel samples from free flowing powders and granules.


The sample volume can be quickly and easily changed by changing removable sample cells.

- 316 stainless steel construction
- Variable sample volumes 5ml to 60ml
- Range of sampler lengths
- Suitable for free flowing powders and granules



Operation

1. Select the correct volume inserts and fit them into the inner rod.
2. Slide the assembled inner rod into the outer tube
3. Insert the sampler into the product, turn the handle to expose the sampling cells to the product.
4. Turn the handle so no sampling inserts are exposed then remove the sampler from the product.
5. Remove the inner rod from the outer sleeve, the individual Sampling Inserts can then be easily removed for analysis.



QUOTE REQUEST ✓

BUY ONLINE ✓

| Stock No. | Length | Diameter | Inserts Required | Material |
|---------------|--------|----------|------------------|---------------------|
| 675 1290K1000 | 1000mm | 40mm | 3 | 316 stainless steel |
| 675 1290K1500 | 1500mm | 40mm | 4 | 316 stainless steel |
| 675 1290K2000 | 2000mm | 40mm | 4 | 316 stainless steel |

Sample Inserts

Inserts are available any volume from 5ml up to 60ml

675 1290W XX Sample Insert where XX = sample volume in ml e.g. 675 1290W25 = insert with 25ml sample volume